

RWD

M1324R


High-Speed Microcentrifuge

A Leading Manufacturer for Pre-clinical Research



M1324R High-Speed Microcentrifuge

....



- ⊗ Maximum capacity: 24 × 1.5/2.0 mL
- ⊗ Maximum speed: 21,130 × g (15,000 rpm)
- ⊗ Temperature control range: -10 °C to +40 °C, only take 8 minutes from room temperature to 4°C (Fast temperature mode); the temperature remains 4°C at the highest speed
- ⊗ Touch screen as well as push button design, for clearer display of the parameters and easier operation
- ⊗ RT mode and Cooling mode, one keyboard toggle
- ⊗ Automatic countdown display mode under unmanned status while spinning
- ⊗ Highly intelligent analysis software: the LED display allows you to preprogram
- ⊗ 12 programs simultaneously for directly use
- ⊗ Multiple abnormal alarm including over-speed, over-temperature, as well as unbalance, ensuring the safety of the operation of the user.
- ⊗ Rear-mounted compressor design, saving space

Technical parameters:

Temperature range	-10 °C to +40 °C
Max. speed	21,130 × g (15,000 rpm)
Max. capacity	24 × 1.5/2.0 mL
Display	4.3-inch touchscreen, 800*480 px
Dimensions (W×D×H)	300mm*514mm*297mm

RWD Life Science Co.,Ltd

Add: Room201, Robeta Building, No.11 High Tech North 4th Road, Nanshan District, Shenzhen City, Guangdong Province, 518057, China
 Tel: (86)-755-86111286 ext 8283 Web: www.rwdstco.com E-mail: sales@rwdstco.com

RWD Life Science Inc.

Add: 6540 Lusk Blvd, Suite C161, San Diego, CA 92121
 Tel: (858)900-5879 E-mail: sales@rwdstco.com

