

*made to measure*


# ACI-01

for perfect chlorinating of silver wires



## Features:

- ⇒ Low-cost, fully automatic system for manufacturing Ag-AgCl electrodes
- ⇒ Makes the fabrication of chlorinated silver wires reliable and easy
- ⇒ Uses an external power supply
- ⇒ Works according to an optimized chlorinating procedure developed at the Max-Planck Institute in Göttingen, Germany
- ⇒ Automatic shut-off
- ⇒ Operation indicated by LED



ACI-01 working scheme

## For more information contact:

**General:**  
**npi electronic GmbH**  
Phone: +49-7141-9730230  
Fax: +49-7141-9730240  
sales@npielectronic.com  
www.npielectronic.com

**North America:**  
**ALA Scientific Instruments**  
Phone: +1-631-393-6401  
Fax: +1-631-393-6407  
sales@alascience.com  
www.alascience.com

**Great Britain:**  
**Scientifica Limited**  
Phone: +44-1825-749933  
Fax: +44-1825-749934  
info@scientifica.uk.com  
www.scientifica.uk.com

**Switzerland:**  
**Science Products Trading AG**  
Phone: +41-43-4880561  
Fax: +41-43-4880562  
info@science-products.com  
www.science-products.ch